

# Impella® Procedure with **SINGLE-ACCESS TECHNIQUE**

## Procedural Steps

Abiomed's Impella CP Introducer\* was recently FDA cleared<sup>1</sup>, and is being released to users. It has a 14 French (F) sheath, which permits placement of the Impella CP's pump (motor and housing). However, after the Impella CP pump is positioned in the ventricle, its smaller Impella 9 F catheter shaft remains in the introducer. The size difference between the 14 F sheath and 9 F catheter shaft permits the insertion of a second introducer sheath (up to 7 F) adjacent to the 9 F catheter shaft, using a new technique. The additional sheath can be used to pass other catheters, including interventional devices. Use of the new technique, called the "single-access approach", results in one access site for both the Impella CP device and an additional catheter.



Only for use with the Impella CP 14 F introducer sheath

# Impella CP® Heart Pump

## SINGLE-ACCESS TECHNIQUE

### Procedural Steps

After you have placed the Impella catheter using access best practices...

#### PROCEDURE STEPS




**Step 1.** Aspirate and flush Impella sheath

**Step 2.** Place needle through diaphragm of Impella 14F sheath at the 10:00 or 2:00 position to allow for placement of 0.035" guidewire

**Step 3.** Secure the Impella catheter to avoid any movement and insert PCI sheath over 0.035" guidewire

**Step 4.** Perform Protected PCI

**Step 5.** Remove PCI sheath at end of procedure, while securing the Impella catheter



Remove Impella and close access site per recommendations

# Single-Access

## TROUBLESHOOTING

with the Impella CP 14 F introducer sheath

### **Bleeding from diaphragm**

- Remove and re-access

### **Bleeding around sheath**

- Tighten Perclose ProGlide® suture

### **Break in peel-away introducer along score lines**

- Remove secondary sheath and peel away the Impella introducer, insert repositioning sheath; gain alternate access for PCI

### **Movement of Impella while advancing sheath**

- Ensure catheter is fixed while advancing sheath

**Abiomed recommends fixating the Impella catheter shaft while advancing sheath**

**Potential Complications:** There are no additional complications related to use of the single-access approach. See the Impella CP Introducer's Instructions for Use found by visiting: <http://www.abiomed.com/impella-device-instructions-for-use> for a complete list of the complications related to its use.

Following current best practice, it is critical to fix the Impella shaft under fluoroscopy while advancing the sheath. Movement of the Impella catheter while advancing the secondary sheath is most common while using non-hydrophilic 7 F sheaths.

Following current best practice, care should be taken when inserting the secondary sheath to prevent sheath damage. Additionally, dilators and catheters should be removed slowly as rapid removal may damage the valve membrane resulting in blood flow through the valve.

When using the Impella 14 F sheath for single-access, a secondary access point is no longer available to manage complications.

**Troubleshooting noted above is the same as current Impella procedures**

## INDICATIONS FOR USE

### High-Risk PCI

The Impella 2.5®, Impella CP® and Impella CP® with SmartAssist® Systems are temporary ( $\leq 6$  hours) ventricular support devices indicated for use during high-risk percutaneous coronary interventions (PCI) performed in elective or urgent, hemodynamically stable patients with severe coronary artery disease, when a heart team, including a cardiac surgeon, has determined high-risk PCI is the appropriate therapeutic option. Use of the Impella 2.5, Impella CP, and Impella CP with SmartAssist Systems in these patients may prevent hemodynamic instability, which can result from repeat episodes of reversible myocardial ischemia that occur during planned temporary coronary occlusions and may reduce peri- and post-procedural adverse events.

### Cardiogenic Shock

The Impella 2.5®, Impella CP®, Impella CP® with SmartAssist®, Impella 5.0®, Impella 5.5® with SmartAssist® and Impella LD® Catheters, in conjunction with the Automated Impella Controller (collectively, "Impella® System Therapy"), are temporary ventricular support devices intended for short term use ( $\leq 4$  days for the Impella 2.5, Impella CP, and the Impella CP with SmartAssist, and  $\leq 14$  days for the Impella 5.0, Impella 5.5 with SmartAssist and Impella LD) and indicated for the treatment of ongoing cardiogenic shock that occurs immediately ( $< 48$  hours) following acute myocardial infarction or open heart surgery or in the setting of cardiomyopathy, including peripartum cardiomyopathy, or myocarditis as a result of isolated left ventricular failure that is not responsive to optimal medical management and conventional treatment measures (including volume loading and use of pressors and inotropes, with or without IABP). The intent of Impella System Therapy is to reduce ventricular work and to provide the circulatory support necessary to allow heart recovery and early assessment of residual myocardial function.

### Important Risk Information for Impella devices

#### CONTRAINDICATIONS

The Impella 2.5, Impella CP, Impella CP with SmartAssist, Impella 5.0, Impella 5.5 with SmartAssist and Impella LD are contraindicated for use with patients experiencing any of the following conditions: Mural thrombus in the left ventricle; Presence of a mechanical aortic valve or heart constrictive device; Aortic valve stenosis/calcification (equivalent to an orifice area of 0.6 cm<sup>2</sup> or less); Moderate to severe aortic insufficiency (echocardiographic assessment graded as  $\geq +2$ ); Severe peripheral arterial disease precluding placement of the Impella System; Significant right heart failure\*; Combined cardiorespiratory failure\*; Presence of an Atrial or Ventricular Septal Defect (including post-infarct VSD)\*; Left ventricular rupture\*; Cardiac tamponade\*

\* This condition is a contraindication for the cardiogenic shock indication only.

#### POTENTIAL ADVERSE EVENTS

Acute renal dysfunction, Aortic valve injury, Bleeding, Cardiogenic shock, Cerebral vascular accident/Stroke, Death, Hemolysis, Limb ischemia, Myocardial infarction, Renal failure, Thrombocytopenia and Vascular injury

In addition to the risks above, there are other **WARNINGS** and **PRECAUTIONS** associated with Impella devices. Visit [www.abiomed.com/important-safety-information](http://www.abiomed.com/important-safety-information) to learn more.

1. See [www.accessdata.fda.gov/SCRIPTS/cdrh/cfdocs/cfPMN/pmn.cfm?ID=K192769](http://www.accessdata.fda.gov/SCRIPTS/cdrh/cfdocs/cfPMN/pmn.cfm?ID=K192769)

\*Catalog numbers: 0052-3056, 0052-3057, 0052-3058

**For further information on single-access with the Impella 14 F sheath, please refer to the Instructions for Use or contact your local team and Abiomed's Clinical Support Center, 1-800-422-8666.**



**ABIOMED, Inc.**  
22 Cherry Hill Drive, Danvers, MA 01923 USA  
Voice: 978-646-1400  
Facsimile: 978-777-8411  
Email: [marketing@abiomed.com](mailto:marketing@abiomed.com)

**Clinical Support**  
**24 hours per day, 7 days a week:**  
1-800-422-8666 (US)