



Leveraging Device-Heart Interaction of the Impella Trans-valvular Pump to Manage Critically Ill Patients

Steven Keller, MD PhD

on behalf of

Elazer Edelman, MD PhD; Brian Chang, PhD; Samantha Polak MD PhD;
Christian Moyer PhD; Qing Tan PhD; Ahmad El Katerji PhD; Noam Josephy
MD

Steven Keller, MD PhD

Research support:

Abiomed (Elazer Edelman, PI)

Scientific and clinical advisory boards:

Abiomed

Paragonix Technologies



Consulting and Other:

CBSET


X-COR Therapeutics

Hemodynamics are critical to managing MCS patients

- Hemodynamic monitoring is critical to managing cardiogenic shock patients^{1,2,3}
- Improved outcomes observed with HD monitoring of in current era of medicine⁴
- PA catheters have inherent limitations⁵ and risks⁶, magnified by infrequent use⁷
- Hemodynamic metrics guide shock protocols^{8,9} but PAC use has plateaued^{1,3,4}



Existing technology carries inherent limitations and risks




- Poor repeatability
- Discreate sampling
- Irrelevant timescale
- Indirect LVEDP (wedge)
- Variable with valve disease
- Risk of adverse events
- Additional access to manage

Marik *Ann of Intensive Care* 2013; Chang, Edelman, et al. *Sci Trans Med* 2018; Chang, Edelman, et al. *Transact on Biomed Eng* 2020

Impella as a Therapeutic Tool → Impella as a Diagnostic Tool



Unique position

1. Indwelling
2. Works in concert with the heart
3. Rotor-motor size and design



Impella LVEDP: Pump responds to variable loading

- Impella pump operates to maintain a fixed rotational speed (e.g. 23k-44k RPM → P1-P8)
- Impella motor current responds to variable load on the device pump
 1. Pressure head
 2. Flow



$$\tau = H * d$$

$H \equiv$ Pressure Head



$d \equiv$ Volumetric Displacement

$$P = \frac{\tau * \omega}{\eta} \quad \eta \equiv \text{Efficiency}$$




$$I_{\text{electrical}} \propto \frac{\tau}{\eta_{\text{combined}}}$$

Impella LVEDP measurement validated in preclinical studies

- Performance evaluated in N=18 acute preclinical studies
- Strong correlation with Millar pigtail reference ($r > 0.80$, bias < 1 mmHg)
- Improved accuracy vs. indirect (wedge) measurements (4 vs. 6 mmHg RMSE)
- Consistent trending in LVEDP over elevated and depressed physiologic states



Impella CO: Novel method to directly measure SVR





D

$$i_{heart} + i_{pump} = \frac{P}{R} + C \frac{dP}{dt}$$
$$P = P_0 e^{-\frac{t}{RC}} + i_{pump} R$$

$$CO = MAP / SVR$$


Impella CO measurement validated in preclinical studies

- Performance evaluated in N=12 acute preclinical studies
- Strong correlation with bolus thermodilution, equivalent measurement accuracy (PE < 30%)
- Improved repeatability vs. individual thermodilution injections (6 vs. 18%)
- Consistent trending in CO over elevated and depressed physiologic states




Currently enrolling in prospective clinical study

SmartPump Study Flow




Data Collection: Synchronized, High Quality HD Reference Data



Design of the Impella *Hemodynamics Platform*

Built on the
Impella CP with **SmartAssist** platform

With **new** (investigational) **features**:




Cardiac Output

Continuous
Accurate
Internal calibration

LVEDP

Continuous
Accurate
Self calibration

Continuous metric trends can improve patient management





Continuous LVEDP,
MAP, and CO

Pump flow vs.
Native Heart CO

High-res,
synchronized
trends


Innovation fueled by academic-industry collaboration



| Criteria No. | Description |
|--------------|---|
| 1 | Successful pacing calibration |
| 2 | The presence of an Impella position alarm |
| 3 | Sustained ventricular arrhythmias requiring cardioversion or temporary pacing |
| 4 | Severe hypotension |
| 5 | Device malfunctions |
| 6 | Device failures |

| Subject No. | 1 | 2 | 3 | 4 | 5 | 6 |
|-------------|-----|-----|-----|-----|-----|-----|
| 1 | 0/0 | 0/0 | 0/0 | 0/0 | 0/0 | 0/0 |
| 2 | 0/0 | 0/0 | 0/0 | 0/0 | 0/0 | 0/0 |
| 3 | 0/0 | 0/0 | 0/0 | 0/0 | 0/0 | 0/0 |
| 4 | 0/0 | 0/0 | 0/0 | 0/0 | 0/0 | 0/0 |
| 5 | 0/0 | 0/0 | 0/0 | 0/0 | 0/0 | 0/0 |
| 6 | 0/0 | 0/0 | 0/0 | 0/0 | 0/0 | 0/0 |
| 7 | 0/0 | 0/0 | 0/0 | 0/0 | 0/0 | 0/0 |
| 8 | 0/0 | 0/0 | 0/0 | 0/0 | 0/0 | 0/0 |
| Total | 0/0 | 0/0 | 0/0 | 0/0 | 0/0 | 0/0 |


*Denominators reflect number of subjects evaluated. Numerators reflect number of subjects that met the safety criteria outlined in Table 2.7
 *0/0 denotes that one (1) pacing attempt at max flow. Pacing was not attempted (due to pump suction)
 *0 denotes minor repositioning after pacing caused by subject's movement, no position alarm was triggered



Mechanisms to measure LVEDP, CO translated to Impella system



**First-in-Human – 8 patients
Evaluated safety and feasibility**

Foundational research in device-heart interactions



Chang et al., STM 2018 Chang et al., IEEE 2019

**Preclinical – 20 validation studies
Proving equivalence to diagnostic devices**



Keller et al., IEEE TBME Keller et al., Annals BME 2020
2020



QUESTIONS

Cited Studies

1. Hernandez et al. Trends in utilization and outcomes of pulmonary artery catheterization in heart failure with and without cardiogenic shock, **JCF** Vol. 25, p.364-371, 2019
2. O'Neill et al. Analysis of Outcomes for 15,259 US patients with acute myocardial infarction cardiogenic shock (AMICS) supported with the Impella device, **AHJ** Vol. 202, p.33-38, 2018
3. O'Neill Adoption of best practices improves survival for patients with AMI-CGS support with Impella, **TCT** 2018, [based on] Abiomed Impella Quality (IQ) Database, US AMI/CGS Apr 2009 – Oct 2017.
4. Ikuta et al. National trends in use and outcomes of pulmonary artery catheters among Medicare beneficiaries, 1999-2013, **JAMA Cardiology** Vol. 2, p.908-913, 2017.
5. Nadeau et al. Misinterpretation of pressure measurements from the pulmonary artery catheter & Limitations of cardiac output measurements by thermodilution, **CASJ** Vol. 33 p.352-363 & p.780-784, 1986.
6. Pandey et al. Use of pulmonary artery catheterization in US patients with heart failure, 2001-2012, **JAMA Internal Medicine** Vol. 176, p.129–132, 2016.
7. Iberti et al. A multi-center study of physicians' knowledge of the pulmonary artery catheter, **JAMA** Vol. 264, p.2928-2932, 1990.
8. Basir et al. Feasibility of early mechanical circulatory support in acute myocardial infarction complicated by cardiogenic shock: The Detroit Cardiogenic Shock Initiative, **CCI** Vol. 91, p.454-461, 2018.
9. Tehrani et al. Standardized team-based care for cardiogenic shock, **JACC** Vol. 73, p.1659-1669, 2019.